

Report 2 - Spectrograph

April 12, 2019

1 Phys 581 Winter 2019

2 Report #2: Live spectrograph data

2.1 Alexander Hickey, 10169582

Note that the contents of this notebook were created and tested in a 64-bit distribution of Windows 10, using Python 3.6.8.

```
In [1]: import sys
        sys.version
```

```
Out[1]: '3.6.8 |Anaconda, Inc.| (default, Feb 21 2019, 18:30:04) [MSC v.1916 64 bit (AMD64)] '
```

```
In [2]: #Import useful libraries
        from IPython.display import Image
```

2.1.1 Introduction

In assignment 6, we explored the construction of a graphical user interface for an Ocean Optics USB2000+ spectrograph. A simple GUI was developed using the Tkinter library in Python, and allowed the user to collect data from the spectrograph with the click of a button, and display this data on screen in a window. Additionally, the user was able to control the integration time of the spectrograph by manually entering the time in a text box. A full discussion on how to interface the spectrograph to a Windows 10 machine can be found in Assignment 6.


While the previous GUI was functional for simple visualization of the data collected by the spectrograph, the user was limited to displaying a single collection at the press of a button. This report presents improvements on the GUI, the most significant one being the implementation of live plotting the data collected by the spectroscope. Additionally, the integration time interface was changed from an entry box to a slider, which defaults to a reasonable setting, and prevents the user from entering an invalid value. Finally, a “max intensity” slider was added, which allows the user to manually adjust the range of the vertical axis on the plot.

2.1.2 The GUI

This section includes several screenshots of the GUI, including the added sliders. These sliders were implemented using the Scale object in tkinter. The following image shows the GUI at startup, before any data has been collected.

```
In [3]: Image(filename="open.JPG",width = 400)
```


Out [3]:



Upon pressing the “Collect” button, the program interfaces with the spectrograph and begins data collection at the desired integration time. The “Collect” button changes to read “Pause”, and upon pressing this button, the plot freezes at the current frame, which can subsequently be manipulated and saved using the matplotlib toolbar. The following image shows the measured spectrum of my computer screen, at the default integration time.

```
In [4]: Image(filename="default.JPG",width = 400)
```


Out [4]:



The following image shows the data collected, again from my computer screen, but for a much larger integration time. In this case, the intensity of the peaks increased greatly, so the vertical axis needed to be adjusted using the “Max intensity” slider.

```
In [5]: Image(filename="increased_integration.JPG",width = 400)
```


Out [5]:



Finally, the next image shows the data, once again collected from my computer screen, but for a much lower integration time.

```
In [6]: Image(filename="decreased_integration.JPG",width = 400)
```

Out[6]:



It is worth noting that the two toolbars work during both the live plotting mode, and paused mode of the application. The full source code for this application can be found in the included SpectrographGUI_Live.py file.

2.1.3 Conclusion

This report presented improvements on the GUI developed in Assignment 6, used to ease the data collection upon interfacing with an Ocean Optics USB2000+ spectrograph. The improvements to this application include the live plotting functionality, as well as the implementation of sliders to control the spectrograph integration time, and the vertical scale of the plot.

With additional time, further functionality could be introduced, such as the ability to save the data file corresponding to the on screen plot while paused. It would also be interesting to spend some time improving the aesthetic of the application, to be more visually appealing to the user.

Appendix: SpectrographGUI_Live.py

```
1  """
2  @author: Alex Hickey
3
4  This program generates a graphical user interface for an Ocean Optics USB2000+
5  spectrograph. The interface includes a collection button, which will interface
6  with the spectrograph and plot the intensity versus wavelength and display
7  the output in real time. There is also a slider to manually set the integration
8  time, as well as the maximum intensity on the plot.
9  """
10
11
12  #Import libraries
13  import tkinter as tk
14  from tkinter import ttk
15  from matplotlib.backends.backend_tkagg import FigureCanvasTkAgg, NavigationToolbar2Tk
16  from matplotlib import pyplot as plt
17  from matplotlib import style
18  import matplotlib.animation as animation
19  import seabreeze.spectrometers as sb
20
21  #Set font, plot style and collection range of spectrograph
22  FONT = ('Verdana', 12)
23  style.use('ggplot')
24  collec_range = (193,896) #Collection range in nanometers!
25  int_range = (0,9500) #Default intensity range
26
27
28
29
30  def get_data(integration_time = 20000):
31      '''
32      This function interfaces with the spectrograph and returns the
33      intensity as a function of wavelength over a given integration time.
34
35      Args:
36          integration_time: Integration time in microsecons
37
38      Return:
39          wlength, intensity: arrays of wavelength and intensity of measurement
40
41      '''
42
43      #Enter loop to continuously update data
44      while True:
45
46          #Connect to spectrometer, set integration time
47          spec = sb.Spectrometer.from_serial_number()
```

```

48         spec.integration_time_micros(integration_time)
49
50         #Retrieve wavelength and intensity arrays
51         wlength = spec.wavelengths()[2:]
52         intensity = spec.intensities()[2:]
53
54         #Close connection
55         spec.close()
56
57         return wlength, intensity
58
59
60
61 class App(tk.Frame):
62     '''
63     Class corresponding to main application. Object is the homepage.
64     '''
65
66     def __init__(self, *args, **kwargs):
67
68         #Initialize homepage
69         tk.Frame.__init__(self, *args, **kwargs)
70         cont = tk.Frame(self)
71         cont.pack(side='top', fill='both')
72
73         #Set application title
74         self.wininfo_toplevel().title('Ocean Optics USB2000+')
75
76         #Attribute to track if animation is running
77         self.running = False
78
79         #Initialize animation attribute
80         self.ani = None
81
82         #Defines integration time label and slider
83         lbl = ttk.Label(cont, text="Integration time ( s ): ", font = FONT)
84         lbl.pack(side=tk.LEFT)
85
86         self.slide = tk.Scale(cont, orient='horizontal', from_= 1000, to= 200000)
87         self.slide.set(20000)
88         self.slide.pack(side=tk.LEFT)
89
90
91         #Defines max intensity label and slider
92         lbl2 = ttk.Label(cont, text=" Max intensity: ", font = FONT)
93         lbl2.pack(side=tk.LEFT)
94
95         self.slide_max = tk.Scale(cont, orient='horizontal', from_= 1000, to_=90000)
96         self.slide_max.set(int_range[1])
97         self.slide_max.pack(side=tk.LEFT)
98
99         #Defines data collection button
100        self.btn = tk.Button(cont, text=' Collect ', command=self.on_click)

```

```

101     self.btn.pack(side=tk.LEFT)
102
103     #Defines the main plot
104     self.fig = plt.Figure(figsize= (10,8))
105     self.ax1 = self.fig.add_subplot(111)
106     self.ax1.set_xlim(*collec_range)
107     self.ax1.set_ylim(*int_range)
108     self.ax1.set_xlabel('Wavelength (nm)',fontsize = FONT[1])
109     self.ax1.set_ylabel('Intensity',fontsize = FONT[1])
110     self.line, = self.ax1.plot([], [], lw=2)
111
112     #Defines the canvas to insert matplotlib plot
113     self.canvas = FigureCanvasTkAgg(self.fig, master=self)
114     self.canvas.draw()
115     self.canvas.get_tk_widget().pack()
116
117     #Inserts the matplotlib toolbar
118     self.toolbar = NavigationToolbar2Tk(self.canvas, self)
119     self.toolbar.update()
120
121
122     def on_click(self):
123         '''
124         Method to decide the action of the collect button
125         '''
126
127         #Start live animation on first click
128         if self.ani is None:
129
130             return self.start()
131
132         #If paused, stop animation
133         elif self.running:
134
135             self.ani.event_source.stop()
136             self.btn.config(text=' Collect ')
137
138         #If resumed, start animation
139         else:
140
141             self.ani.event_source.start()
142             self.btn.config(text=' Pause ')
143
144         #Change status
145         self.running = not self.running
146
147     def start(self):
148         '''
149         Method to start/resume animation
150         '''
151
152         #Signal that animation is running
153         self.running = True

```



```

154
155     #Change button label
156     self.btn.config(text=' Pause ')
157
158     #Start animation
159     self.ani = animation.FuncAnimation(self.fig,
160                                       self.update_graph,
161                                       interval=int(self.slide.get())/1000,
162                                       repeat=False)
163     self.ani._start()
164
165
166     def update_graph(self, i):
167         """
168         Method to update the spectrograph plot with current data.
169         """
170
171         #Retrieve updated integration time
172         int_time = int(self.slide.get())
173
174         #Update plot
175         self.line.set_data(*get_data(int_time))
176         self.ax1.set_ylim(0,self.slide_max.get())
177
178         return self.line,
179
180
181
182     def main():
183         """
184         Compile and run main application.
185         """
186
187         #Create tk object
188         root = tk.Tk()
189
190         #Set application icon
191         root.iconbitmap(default = 'light_icon.ico')
192
193         #Compile and execute application
194         App(root).pack()
195         root.mainloop()
196
197
198     if __name__ == '__main__':
199         main()

```